

## Protocol report

NucleoMag\_DNA\_RNA\_Water\_Flex\_Rev02

5/20/2021 5:10:03 PM-04:00

1/5

# General info

## Protocol information

Protocol name	NucleoMag_DNA_RNA_Water_Flex_Rev02
Modified by	admin
Kit name	NucleoMag DNA/RNA Water
Description	Isolation of DNA and RNA from water.
	NucleoMag DNA/RNA Water REF: 744220.1 /.4
	Rev. 1811/002

## Protocol report

NucleoMag\_DNA\_RNA\_Water\_Flex\_Rev02

5/20/2021 5:10:03 PM-04:00

2/5

## Sample layout

## Reagent info

Tipcomb		96 standard plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
-	-	-	-

Binding		96 DW plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
Cleared lysate	450	-	Reagent
NucleoMag B-Beads	25	-	Reagent
Binding Buffer MWA2	475	-	Reagent

1st Wash MWA3		96 DW plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
Wash Buffer MWA3	850	-	Reagent

2nd Wash MWA3		96 DW plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
Wash Buffer MWA3	850	-	Reagent


3rd Wash MWA4		96 DW plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
Wash Buffer MWA4	850	-	Reagent

Elution		96 standard plate	
Name	Well volume [µl]	Total reagent volume [µl]	Type
RNase-free H2O	100	-	Reagent

## Dispensed reagents

The protocol does not contain dispensed reagents

## Steps data




	Tip1	96 DW tip comb	
	Pick-Up	Tipcomb	
	Binding	Binding	
	Beginning of step	Precollect	No
		Release beads	Yes
	Mixing / heating:	Shake 1 time, speed	00:00:30, Medium
		Shake 2 time, speed	00:00:30, Bottom mix
		Loop count	10
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	3
	1st Wash	1st Wash MWA3	
	Beginning of step	Precollect	No
		Release beads	Yes
	Mixing / heating:	Shake 1 time, speed	00:00:30, Medium
		Shake 2 time, speed	00:00:30, Bottom mix
		Loop count	2
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	3
	2nd Wash	2nd Wash MWA3	
	Beginning of step	Precollect	No
		Release beads	Yes
	Mixing / heating:	Shake 1 time, speed	00:00:30, Medium
		Shake 2 time, speed	00:00:30, Bottom mix
		Loop count	2
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	3
	3rd Wash	3rd Wash MWA4	
	Beginning of step	Precollect	No
		Release beads	Yes
	Mixing / heating:	Shake 1 time, speed	00:00:30, Medium
		Shake 2 time, speed	00:00:30, Bottom mix
		Loop count	2
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	3

## Protocol report

NucleoMag\_DNA\_RNA\_Water\_Flex\_Rev02

5/20/2021 5:10:03 PM-04:00

5/5

	Dry	3rd Wash MWA4	
		Dry time	00:10:00
		Tip position	Outside well / tube
	Elution	Elution	
		Beginning of step	Precollect Release beads No Yes
	Mixing / heating:	Shake 1 time, speed	00:00:30, Medium
		Shake 2 time, speed	00:00:30, Bottom mix
		Loop count	5
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	5
		Collect time [s]	25
	Leave	Tipcomb	

## Lot info

No lot numbers have been defined.