



## Reagent info



Sample Plate		96 DW plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
Sample Lysate	200	-	Sample
ZymoBIOMICS MagBinding Beads	50	-	Reagent
ZymoBIOMICS MagBinding Buffer	600	-	Reagent
ZymoBIOMICS W1		96 DW plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
ZymoBIOMICS Wash 1	500	-	Reagent
ZymoBIOMICS W2 1		96 DW plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
ZymoBIOMICS Wash 2	900	-	Reagent
ZymoBIOMICS W2 2		96 DW plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
ZymoBIOMICS Wash 2	900	-	Reagent
Elution Plate		96 DW plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
ZymoBIOMICS DNAse/RNAse Free Water	50	-	Reagent
Tip Plate		96 standard plate	
Name	Well volume [μl]	Total reagent volume [μl]	Type
-	-	-	-

## Dispensed reagents

The protocol does not contain dispensed reagents

## Steps data

	Tip	96 DW tip comb	
	Pick-Up	Tip Plate	
	Binding	Sample Plate	
	Beginning of step	Precollect	No
		Release time, speed	00:00:30, Fast
	Mixing / heating:	Mixing time, speed	00:10:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	2
		Collect time [s]	5
	MagWash 1	ZymoBIOMICS W1	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:05:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	2
		Collect time [s]	5
	MagWash 2-1	ZymoBIOMICS W2 1	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:05:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	2
		Collect time [s]	5
	MagWash 2-2	ZymoBIOMICS W2 2	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Fast
	Mixing / heating:	Mixing time, speed	00:05:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	2
		Collect time [s]	5
	Dry	ZymoBIOMICS W2 2	
		Dry time	00:10:00
		Tip position	Outside well / tube

	Elution	Elution Plate	
	Beginning of step	Precollect	No
		Release time, speed	00:00:10, Slow
	Mixing / heating:	Mixing time, speed	00:10:00, Bottom mix
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	2
		Collect time [s]	5
	Leave	Tip Plate	

## Lot info

No lot numbers have been defined.